


Ordering options:

- a) SCI-FERM-COMP/HSV2
- b) SCI-FERM-COMP/ HS-WP1T-20L



Features:

- Oil less compressor/vacuum pump
 - Unique Teflon coated Piston/diaphragm
 - Continuously rated, rugged 1f motor
 - Permanently lubricated sealed ball bearings
 - Practically maintenance free
 - Ultra light die cast aluminum components
 - Minimal weight, low noise
 - Balanced eccentrics for smooth low vibration operation
- Built in micro-air filter



Technical Specifications

Model	HSV-2	HS-WP1T-20L		
Type	Diaphragm vaccum	Comp with 20L tank		
Motor HP	1/3	3/4		
Volts/Hertz/Phase	230/50/1	230/50/1		
Weight	9.0 Kgs	40.0 Kgs		
Dimensions	9" X 4.5" X 13"	18" X 9" X 20"		
Compressor performance	PSIG	LPM	PSIG	LPM
	0	70	0	54
	10	60	10	52
	20	45	20	47
	30	26	30	42
	40	10	40	38
			50	34
			60	32
			80	28
			100	22
Vacuum Performance	IN mmHg	LPM	IN mmHg	LPM
	0	50	-	-
	5	39	-	-
	10	28	-	-
	15	17	-	-
	20	12	-	-
	25	5	-	-