

Pressure pad Analytical pump

Product Summary

Our first generation pump based on stepper motor drive with logic sequencer.

These pumps offer extremely high pressure & many different tube compositions can be used.

- Digital Logic sequencer based or
- Stepper motor drive for uniform rpm
- Pump speed doesn't vary with power fluctuations
- Eight roller, ultra low pulsation
- Independent tube pressure control
- 0-100 rpm in steps of 10 (or custom flow)
- 1 or 2 channel pressure pad
- Molybdenum disulphide Nylon rollers



High pressure analytical pressure pad peristaltic pump

Ordering options:

Part number

SCI/PP-ANAL/SM-LS/8X

SCI/PP-ANAL/SM-MICR/8X

Description

Analytical 8 Roller 1 or 2 channel pressure pad pump system, Logic sequencer based stepper motor drive

Analytical 8 Roller 1 or 2 channel pressure pad pump system, microcontroller based, Rs232 interface, DS1307 RTC, 16X2 Alphanumeric LCD with 4 button keypad, built in timer & external trigger