

Ordering options:

- a) SCI-FERM-S.I.P/AUTO-CLR-TCH5.7/5,9,18L
- b) SCI-FERM-S.I.P/SEMI-AUTO/5,9,18L

Model	5L	9L	18L
Type	Auto feed	Auto feed	Auto feed
Water level	3.5L approx	3.6L approx	6.3L
Vessel capacity	5.4L approx	9.0L approx	18.0L approx
Power	240V-2000W	240V-5400W	240/440V-9000W
Operating pressure	3.0-3.25 bars	3.0-3.25 bars	3.0-3.25 bars
Construction	Stainless steel	Stainless steel	Stainless steel
Control: Micro-controller	ATMEGA16/32	ATMEGA16/32	ATMEGA16/32
Pressure sensors: 2 Nos	Honeywell-Sensym	Honeywell-Sensym	Honeywell-Sensym
Solenoids	3 steam, 2 Air	3 steam, 2 Air	3 steam, 2 Air

Note: Connections are via wire-braided hoses, SS316 compression fittings & touch screen only available on export systems

